

## Operator Preventive Maintenance Checks and Services

4110-01-287-7111

Refrigerator, Solid State, Biological, Model, RCB42P

[B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

ITEM NO	INTERVAL	ITEM TO BE INSPECTED AND PROCEDURE	IS NOT MISSION CAPABLE IF:
1	B, A	<b>Refrigerator, Solid State, Biological</b> a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.	Missing components preclude operation of the refrigerator.
	B	b. Inspect the unit for dust, dirt, or damage. Inspect power cord for cracks or tears.	Damage or deteriorated components prevent the operation of the unit, the power cord is cracked or frayed, or wires are not covered by the cord insulation.
	B	c. Assemble the refrigerator power interconnections as directed in the manufacturer's literature. Ensure the proper power source is correct.	The proper power connectors cannot be used.
	B	d. Verify that the hinges and catches are tightly fixed.	The hinges and/or catches are not functional.
	B	e. Verify that the lid seals.	The lid does not seal.
	B	f. Clean the refrigerator as directed in the manufacturer's literature.	
2	B	<b>DC – AC Operation</b> Conduct the operating procedures as directed by the operator's manual.	The refrigerator cannot be started.