

## Repairer Preventive Maintenance Checks and Services

4110-01-287-7111

Refrigerator, Solid State, Biological, Model, RCB42P

[M-Monthly, Q-Quarterly, S-Semiannually, and A-Annually]

ITEM NO	INTERVAL	ITEM TO BE INSPECTED AND PROCEDURE	IS NOT MISSION CAPABLE IF:
1	S	<b>Refrigerator, Solid State, Biological</b> Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand.	Missing components or accessories prevent the operation of the refrigerator.
2	S	<b>Maintaining Refrigerator</b> a. Verify that the refrigerator functions as directed by the Operators Preventative Maintenance Checks and Services manual.  b. Verify that the unit functions on AC current.  c. Verify that the unit functions on DC current.  d. Check screw connections as directed by the manufacturer's service manual.  e. If necessary, conduct the "ACCU" as directed by the manufacturer's service manual.  f. Verify the electrical safety.	The refrigerator cannot function on AC.  The refrigerator cannot function on DC.          The refrigerator fails any of the electrical safety tests.