

Operator Preventive Maintenance Checks and Services

6630-01-524-4220

Mini VIDAS

B-Before Operation, D-During Operation, A-After Operation, Q-Quarterly, and S-Semiannually]

ITEM NO	INTERVAL	ITEM TO BE INSPECTED AND PROCEDURE	IS NOT MISSION CAPABLE IF:
1	B, S	Mini VIDAS a. Conduct an inventory to ensure that the items listed on the Equipment Parts and Accessories List are on hand. b. Check the electrical power cord for cuts, fraying, or deterioration. c. Perform "Start-up" procedures in accordance with operating instructions.	Missing parts or accessories preclude operation of the refrigerator. The power cord is cracked or frayed, wires are not covered by the cord insulation, or the damage prevents the mini VIDAS from working. Mini VIDAS fails to start up.
2	B, S	Display Screen	Display fails.
3	B, A, S	Keypad	Unable to select modes.
4	B, S	Internal Printer	Will not print.
5	B, D, A	Equipment Care Ensure "Equipment Care" is conducted as directed by the manufacturer's literature.	