

## (continued) Appendix A. Operator PMCS

Table A-2. TM 8-4110-001-24&amp;P, Refrigerator, Mechanical Blood Bank

ITEM NO	INTERVAL					ITEM TO BE INSPECTED AND PROCEDURE	EQUIPMENT IS NOT READY/AVAILABLE IF:
	B	D	A	Q	S		
1	X				X	<b>Refrigerator.</b> a. Check for broken, worn, or damaged switches, indicators, and displays on the control panel. b. Check the electrical power cord for cuts, fraying, or deterioration. c. Perform a complete operational checkout of the refrigerator. Ensure that the 9-volt standby battery (located on top of the monitor module rear protective/dust cover) is operational.	Damage prevents proper operation of the refrigerator.  Worn, damaged, or deteriorated power cord prevents safe operation of the refrigerator.  Monitor does not function as specified in the operating procedures.
2	X		X		X	<b>Doors.</b> a. Verify that the doors close and seal properly. b. Inspect the door hinges for loose or missing hardware.	Defective door gasket prevents proper operation of the refrigerator.  Loose or missing hardware prevents proper operation or causes a safety hazard.
3	X				X	<b>Drawers.</b> Ensure that the drawers are unobstructed and move freely. Inspect for dirt and other foreign substances.	Obstructed or damaged drawers prevent refrigerator doors from sealing properly.
4					X	<b>Condensing Unit.</b> Inspect fan condensing unit for damage and dust.	Damage prevents proper operation.