

Table A-20. TM 8-6530-005-24&amp;P, Ultrasonic Cleaner, Mobile

ITEM NO	INTERVAL					ITEM TO BE INSPECTED AND PROCEDURE	EQUIPMENT IS NOT READY/AVAILABLE IF:
	B	D	A	Q	S		
1	X	X	X			<b>Cabinet (exterior).</b> a. Check for loose or missing hardware. Check for chemical deposits. b. Check louvers for an accumulation of dust.	Loose or missing hardware prevents operation or creates a safety hazard. Unit overheats or does not operate correctly.
2	X					<b>Fuseholder.</b> Disconnect the unit from electrical power. Check the fuseholder for tightness.	Fuseholder is loose or damaged.
3	X	X				<b>Timer (control).</b> Check for smooth rotary movement. Check for counterclockwise movement during operation.	Damaged timer prevents operation of the unit.
4	X	X				<b>Lamp.</b> Verify that the lamp illuminates during operation.	Lamp does not illuminate.
5	X					<b>Electrical cable.</b> Inspect the electrical cable and connector for damage or deterioration.	Damage prevents operation or causes a safety hazard.
6	X	X				<b>Drain valve.</b> Verify that the drain valve opens and closes properly.	Valve does not open and close properly.