

(continued) Appendix A. Operator PMCS

Table A-21. TM 8-6530-007-24&amp;P, Solution Warming Cabinet, Mdl 5520

ITEM NO	INTERVAL					ITEM TO BE INSPECTED AND PROCEDURE	EQUIPMENT IS NOT READY/AVAILABLE IF:
	B	D	A	Q	S		
1	X					<b>Warming cabinet.</b> a. Visually inspect for damage to the panels, controls, indicators and door. b. Inspect door gasket for dirt or loose fastening screws.	The unit does not operate properly or a safety hazard exists. Damage or dirt prevents door from sealing properly.
2	X					<b>Control module.</b> Verify its performance by following the initial start-up procedures and the functional testing procedures.	The unit does not control the temperature.
3	X					<b>Electrical power cable.</b> Inspect for damage, deterioration, or fraying.	The cable is damaged or an electrical hazard exists.
4	X		X			<b>Shelf filter.</b> Turn power off and inspect filter.	The filter is damaged or clogged.
5	X		X			<b>Air vents.</b> Inspect the air circulating vents for an excessive accumulation of dust.	The collection of dust may result in a fire.