

Table A-8. TM 8-6515-006-24&P, Endoscopic Light, Mdl 52-1201

ITEM NO	INTERVAL					ITEM TO BE INSPECTED AND PROCEDURE	EQUIPMENT IS NOT READY/AVAILABLE IF:
	B	D	A	Q	S		
1	X		X			Power cable assembly. Check for a frayed electrical cable or a damaged connector.	Damaged or broken power cable assembly prevents operation.
2	X	X	X			Endoscopic instrument light. <i>a.</i> Ensure that the multiport faceplates rotate smoothly in either direction. <i>b.</i> Intensity controls do not turn the lamps off or do not vary the intensity of the illumination. <i>c.</i> Verify that at least one "REPLACE SPARE LAMP" indicator is lit. <i>d.</i> Ensure that an airflow can be detected from the rear exhaust fan vents.	Defective operating controls prevent operation with specific fiberoptic cable ports. Defective operating controls prevent proper operation. More than one spare lamp cartridge assembly is unserviceable. The cooling fans are obstructed or inoperable.