

## (continued) Appendix B. Repairer PMCS

Table B-18. TM 8-6520-004-14&amp;P, Dental Operating Chair and Stool Unit, Portable, Mdl CM-185

ITEM NO	INTERVAL					ITEM TO BE INSPECTED AND PROCEDURE	EQUIPMENT IS NOT READY/AVAILABLE IF:
	B	D	A	Q	S		
1					X	<b>Dental chair.</b> <i>a.</i> Verify that components have been inventoried by operator/user personnel. <i>b.</i> Verify proper assembly and operation of the dental chair. <i>c.</i> Check for any oil on the base plate, around the base covers, or in the container bottom section. <i>d.</i> Inspect the upholstery snaps for proper fastening. <i>e.</i> Check operation of the hydrolok assembly in the chair back.	Missing components prevent operation of the dental chair.  Incorrect assembly or operation of the dental chair prevents the treatment of patients.  The loss of oil in the hydraulic system prevents operation of the dental chair.  Damaged or missing upholstery snaps prevent use of the dental chair.  An inoperable hydrolok assembly prevents use of the dental chair.
2					X	<b>Stool.</b> Verify proper assembly and operation of the stool.	Incorrect assembly or operation of the stool prevents the treatment of patients.
3					X	<b>Case.</b> <i>a.</i> Inspect for cracks, major dents, or puncture holes. <i>b.</i> Check for bent, loose, or broken latches.	Case damage prevents protective storage or safe movement of the unit.  Unserviceable latches prevent safe movement.