

SigmaPace™ 1000 External Pacemaker Analyzer



Overview

The SigmaPace™ 1000 is a powerful handheld analyzer with a comprehensive range of test suites, measurement algorithms and test loads – as specified by the external pacemaker manufacturers – to fulfill your testing requirements both quickly and efficiently.

Standard Features

- Unique full-featured biomedical test product
- Tests for both transcutaneous and transvenous external pacemakers
- Full range of user-selectable measurement algorithms and test loads for external pacemakers
- Dual-channel signal acquisition for capturing synchronous AV-sequential pulse data
- Interactive pacemaker and ECG simulation with 5-lead output
- Large 21-character x 8-line alphanumeric LCD readout
- Exclusive readout "HOLD" function
- Exclusive SigmaPace™ 1000 test features: DC static/dynamic leakage, and current drain

[Biomedical](#) > [Performance Analyzers](#) > [External Pacemaker](#) >

SigmaPace™ 1000 External Pacemaker Analyzer

Specifications	
Modes of Operation	Manual Remote (via standard RS232 serial port)
User Interface	Display: - LCD readout - 21 characters X 8 lines - Brightness/viewing angle adjustment Keys: Eight push buttons (F-2, F-3 [UP arrow], F-4 [UP arrow], two DOWN arrows, ESCAPE, and ENTER)
Transcutaneous Pacer Tests:	Output Pulse Measurement Current: - Ranges: 4.00 mA to 9.99 mA; 10.0 mA to 99.9 mA; 100 mA to 250 mA - Accuracy: $\pm 2\%$ of reading or $\pm 50\ \mu\text{A}$ (whichever is greater) Rate: - Ranges: 5.0 PPM to 99.9 PPM; 100 PPM to 300 PPM - Accuracy: $\pm 0.5\%$ of reading or $\pm 0.3\ \text{PPM}$ (whichever is greater) Width: - Ranges: 1.00 mS to 9.99 mS; 10.0 mS to 99.9 mS - Accuracy: $\pm 0.5\%$ of reading or $\pm 14\ \mu\text{S}$ (whichever is greater) Energy: - Ranges: 1 μJ to 999 μJ ; 1 mJ to 999 mJ; 1.00 J to 1.99 J - Accuracy: 5 % of reading/computation Demand and Asynchronous Mode Tests

	Waveform (Physiological Simulation): Amplitude: Modes of Operation: Active Outputs: Pacemaker Compatibility: Amplitude Sensitivity Test Selections: R, S and T Waves Waveforms: Amplitude: Width: 5-mS steps (20 to 95 mS) Active Outputs: Pacemaker Compatibility: Noise Immunity/Line Frequency Test Waveform: Amplitude Testload Output: 5-lead ECG Output: Active Outputs: Paced Refractory Period Test (PRP) Range: Physiological Simulation: Outputs: Pacemaker Compatibility: Sensed Refractory Period Test (SRP) Range: Accuracy: Physiological Simulation: Active Outputs: Pacemaker Compatibility:	<ul style="list-style-type: none"> - Normal sinus rhythm (NSR) - Complete P-QRS-T complex - 1.0mVpeak (lead I) - Underdrive: NSR @ 85 % of pulse interval/rate - Overdrive: NSR @ 115 % of pulse interval/rate - Auxiliary Control: NSR adjustable in 1-BPM increments - Auxiliary Rate Range: Underdrive 10 BPM (min); overdrive 300 BPM (max) - 5-lead ECG; ventricular test load; high-level ECG jack - Pulse Rates: 30 to 200 PPM - Intended Types: Demand: VVI (pace and sense); async: VOO (pace) - Rate: 30 to 200 PPM - Test Loads: (30) 50 Ω to 1550 Ω in 50-Ω steps - Square (SQU); triangle (TRI); haversine (HSN); sine square (SSQ) - Range: 0.05 mVpeak to 5.0 mVpeak - Accuracy: ± 5 % of setting - Resolution: 0.05-mV steps (0.05 to 0.95 mVpeak); 0.50-mV steps (1 to 5 mVpeak) - Range: 0.15 mS to 300 mS - Accuracy: ± 5 % of setting - Selections: 50 - Resolution: 0.05-mS steps (0.15 to 0.95 mS); 1-mS steps (1 to 19 mS); - 25-mS steps (100 to 300 mS) - 5-lead ECG; ventricular test load; high-level ECG jack - Pulse Rates: 30 to 200 PPM - Intended Type: VVI (pace and sense) - Sine wave - Frequency: 50 and 60 Hz - Accuracy: 0.5 Hz - Range: 0 (OFF) to 10 mVpeak-to-peak - Accuracy: ± 5 % of setting - Resolution: 0.5-mVpeak-to-peak steps - Settings: (30) 50 Ω to 1550 $\Omega \pm 1$ % - Range: 0 (OFF) to 10 mVpeak-to-peak - Accuracy: ± 5 % of setting - Resolution: 0.50-mV steps - Reference: Lead I (RA to LA) - 5-lead ECG; ventricular test load - 20 mS to 500 mS - Accuracy: 5 % of reading or 1 mS (whichever is greater) - Selection: Single pulse, R Wave, SSQ - Pulse Width: 40 mS - 5-lead ECG; ventricular test load - Pulse Rates: 30 to 200 BPM - Intended Type: VVI (pace and sense) - 15 mS to 500 mS - ± 5 % of reading or ± 1 mS (whichever is greater) - Selection: Double pulse, R Wave, SSQ - Pulse Width: 40 mS - Amplitude: 1 mVpeak lead I - 5-lead ECG; ventricular test load - Pulse Rates: 30 to 200 BPM - Intended Type: VVI (pace and sense)
--	---	---

	<p>Test Loads</p> <p>Transcutaneous Pacer:</p> <ul style="list-style-type: none"> - Selections: (31) 50 Ω to 1550 Ω in 50-Ω steps - Accuracy: $\pm 1\%$ of selection - Power Rating: 5 W (average); 40 W (peak) @ 1000 Ω <p>Input Defibrillation Protection:</p> <ul style="list-style-type: none"> - Type: Internal spark gap - Episode Limit: 5 pulses @ 360 J (10 seconds min between discharges) - Life Limit: 250 pulses @ 360 J
Transvenous Pacer Tests:	<p>Output Pulse Measurement</p> <p>Current:</p> <ul style="list-style-type: none"> - Ranges: 0.05 mA to .999 mA (available single channel only); 1.00 mA to 9.99 mA; 10 mA to 30 mA - Accuracy: $\pm 2\%$ of reading or $\pm 50\ \mu\text{A}$ (whichever is greater) - Polarity Indicator: + or - <p>Rate:</p> <ul style="list-style-type: none"> - Ranges: 10.0 PPM to 99.9 PPM; 100 PPM to 999 PPM - Accuracy: $\pm 0.5\%$ or 0.3 PPM (whichever is greater) <p>Width:</p> <ul style="list-style-type: none"> - Ranges: 0.020 mS to .999 mS; 1.00 mS to 9.99 mS; 10.0 mS to 99.9 mS - Accuracy: 0.5% or $\pm 14\ \mu\text{S}$ (whichever is greater) - Resolution: $\pm 1\ \text{LSD}$ or $\pm 4\ \mu\text{S}$ (whichever is greater) <p>Voltage:</p> <ul style="list-style-type: none"> - Ranges: (available single channel only) 0.050 Vpeak to .999 Vpeak; 1.00 Vpeak to 9.99 Vpeak; 10 Vpeak to 30 Vpeak - Accuracy: $\pm 2\%$ of reading or $\pm 0.05\ \text{Vpeak}$ (whichever is greater) - Polarity Indicator: + or - <p>Energy:</p> <ul style="list-style-type: none"> - Ranges: (available single channel only) 1 nJ to 999 nJ; 1 μJ to 999 μJ - Accuracy: $\pm 5\%$ of reading/computation <p>Display Formats:</p> <p>Atrial channel only; ventricular channel only; both A + V channels</p> <p>AV Interval (Delay Time) Measurement Ranges:</p> <p>10.0 mS to 99.9 mS; 100 mS to 999 mS</p> <p>Start Point:</p> <p>Atrial pulse leading edge</p> <p>Stop Point:</p> <p>Ventricular pulse leading edge</p> <p>Accuracy:</p> <p>1 % of reading/computation</p> <p>Demand/Async Mode Tests</p> <p>Channels:</p> <p>Single and dual</p> <p>Waveform:</p> <p>Sine square (SSQ)</p> <p>Atrial Output:</p> <p>Simulated P Wave</p> <ul style="list-style-type: none"> - Width: 30 mS - Amplitude: 2.0 mVpeak <p>Vent Output:</p> <p>Simulated R Wave</p> <ul style="list-style-type: none"> - Width: 40 mS - Amplitude: 2.5 mVpeak <p>AV Interval:</p> <p>90 mS (fixed)</p> <p>Interactive Simulated Rates:</p> <ul style="list-style-type: none"> - Default Settings: Underdrive = NSR @ 85 % of pulse interval/rate; overdrive = NSR @ 115 % of pulse interval/rate - Manual: NSR simulations adjustable in 1-BPM increments - Limits: Underdrive (min) = 10 BPM; overdrive (max) = 300 BPM <p>Output:</p> <p>Ventricular channel test load; atrial channel test load</p> <p>Pulse Rate:</p> <p>30 to 200 PPM</p> <p>Intended Pacemaker Types:</p> <ul style="list-style-type: none"> - Demand: VVI (V-channel pace and sense); AAI (A-channel pace and sense); - Async/Continuous: VOO (V-channel pace and sense); AOO (A-channel pace and sense); DOO (dual-channel pace and sense) <p>DDD (dual-channel pace and sense)</p> <p>Amplitude Sensitivity Test Operation:</p> <p>Single-channel operation only (atrial or ventricular)</p> <p>Atrial Channel (Physiological)</p> <ul style="list-style-type: none"> - Selection: P Wave

	<p>Simulation):</p> <ul style="list-style-type: none"> - Rate: 30 to 120 BPM - Timing: Waveform delayed by 80 % of the pulse-to-pulse interval or 400 mS (whichever is shorter) - Active Output: Atrial test load <p>Available Test Loads: Waveform Selections:</p> <p>200 Ω, 500 Ω (default setting) and 1000 $\Omega \pm 1\%$ Square (SQU); triangle (TRI); haversine (HSN); sine square (SSQ) (default setting); asymmetrical triangle (ISO) – fixed width: 2 mS rise time/13 mS fall time</p> <p>Sensitivity Waveform Amplitude:</p> <p>- Test Load Selection:</p> <p>500 Ω (default setting) -- Range: 0.05 mVpeak to 50.0 mVpeak -- Accuracy: $\pm 5\%$ of setting -- Resolution: 0.05 mVpeak (0.05 to 0.95 mVpeak); 0.50 mVpeak (1 to 50 mVpeak)</p> <p>- Test Load Selection:</p> <p>200 Ω -- Range: 0.05 mVpeak to 20.0 mVpeak -- Accuracy: $\pm 5\%$ of setting -- Resolution: 0.05 mVpeak (0.05 to 0.95 mVpeak); 0.50 mVpeak (1 to 20 mVpeak)</p> <p>- Test Load Selection:</p> <p>1000 Ω -- Range: 0.05 mVpeak to 100 mVpeak-to-peak -- Accuracy: $\pm 5\%$ of setting -- Resolution: 0.05 mVpeak (0.05 to 0.95 mVpeak); 0.50 mVpeak (1.0 to 49.5 mVpeak); 05.0 mVpeak (50 to 100 mVpeak) - Default Setting: 2.0 mVpeak - Range: 0.15 mS to 95.0 mS - Accuracy: $\pm 5\%$ of setting - Resolution: 0.05 mS (0.15 mS to 0.95 mS); 1 mS (1 mS to 19 mS); 5 mS (20 mS to 95 mS)</p> <p>Widths:</p> <p>Intended Pacemaker Types:</p> <p>AAI (atrial pace and sense) - Pulse Rates: 30 to 200 PPM - Selection: R Wave, S Wave and T Wave</p> <p>Ventricular Channel (Physiological Simulation):</p> <ul style="list-style-type: none"> - Rate: 30 to 120 BPM - Timing: Waveform delayed from the ventricular demand pacemaker pulse by 80 % of the pulse-to-pulse interval or 400 mS (whichever is shorter) - Active Output: Selected ventricular test load <p>Waveform Selections:</p> <p>Square (SQU); triangle (TRI); haversine (HSN); sine square (SSQ) (default setting); asymmetrical triangle (ISO) – fixed width: 2 mS rise time/13 mS fall time</p> <p>Available Test Load(s): Amplitude:</p> <p>- Pacer Load Selection:</p> <p>500 Ω -- Range: 0.05 to 50.0 mVpeak -- Accuracy: $\pm 5\%$ of setting -- Resolution: 0.05 mVpeak (0.05 to 0.95 mVpeak); 0.50 mVpeak (1 to 50 mVpeak)</p> <p>- Pacer Load Selection:</p> <p>200 Ω -- Range: 0.05 mVpeak to 20.0 mVpeak -- Accuracy: $\pm 5\%$ of setting -- Resolution: 0.05 mVpeak (0.05 to 0.95 mVpeak); 0.50 mVpeak (1 to 20 mVpeak)</p> <p>- Pacer Load Selection:</p> <p>1000 Ω -- Range: 0.05 mVpeak to 100 mVpeak-to-peak -- Accuracy: $\pm 5\%$ of setting -- Resolution: 0.05 mVpeak (0.05 to 0.95 mVpeak); 0.50 mVpeak (1.0 to 49.5 mVpeak); 05.0 mVpeak (50 to 100 mVpeak) - Default Setting: 2.5 mVpeak - Range: 0.15 mS to 300 mS - Accuracy: $\pm 5\%$ of setting</p> <p>Widths:</p>
--	---

	<p>Intended Pacemaker Type(s):</p> <p>Noise Immunity Test</p> <p>Channels:</p> <p>Waveform:</p> <p>Frequency:</p> <p>Accuracy:</p> <p>Active Output(s):</p> <p>Output Selections:</p> <p>ECG Signal:</p> <p>Atrial:</p> <p>Ventricular:</p> <p>Noise Line Amplitude:</p> <p>- Pacer Load Selection:</p> <p>- Pacer Load Selection:</p> <p>- Pacer Load Selection:</p> <p>Refractory Period Test (Atrial Channel)</p> <p>Test Selections:</p> <p>Period:</p> <p>Accuracy:</p> <p>Resolution:</p> <p>Physiological Simulation:</p> <p>Additional Waveform Selections:</p> <p>Amplitude:</p> <p>Width:</p> <p>Active Outputs:</p> <p>Intended Pacemaker Types:</p> <p>Pacemaker Rates:</p> <p>Available Test Load:</p> <p>Refractory Period Test (Ventricular Channel)</p> <p>Test Selections:</p> <p>Period:</p> <p>Accuracy:</p> <p>Resolution:</p> <p>Display Format:</p> <p>Physiological Simulation:</p>	<p>- Resolution: 0.05 mS (0.15 to 0.95 mS); 1 mS (1 to 19 mS); 5 mS (20 mS to 95 mS); 25 mS (100 mS to 300 mS)</p> <p>VVI (atrial pace and sense only)</p> <p>- Pulse Rates: 30 to 200 PPM</p> <p>Single, atrial or ventricular only</p> <p>Sine wave</p> <p>50 and 60 Hz</p> <p>± 0.5 Hz</p> <p>Atrial- or ventricular-channel test load</p> <p>Atrial channel only; ventricular channel only</p> <p>ECG signal can be added to the selected channel.</p> <p>1 mVpeak</p> <p>1.5 mVpeak</p> <p>500 Ω</p> <p>-- Range: 0 (OFF) to 100 mVpeak-to-peak</p> <p>-- Accuracy: ± 5 % of setting</p> <p>-- Resolution: 5-mVpeak-to-peak steps</p> <p>200 Ω</p> <p>-- Range: 0 (OFF) to 40 mVpeak-to-peak</p> <p>-- Accuracy: ± 5 % of setting</p> <p>-- Resolution: 5-mVpeak-to-peak steps</p> <p>1000 Ω</p> <p>-- Range: 0 (OFF) to 200 mVpeak-to-peak</p> <p>-- Accuracy: ± 5 % of setting</p> <p>-- Resolution: 5-mVpeak-to-peak steps</p> <p>Paced refractory period; sensed refractory period</p> <p>20 to 500 mS</p> <p>± 5 % of reading (or ± 1 mS, whichever is greater)</p> <p>± 1 LSD</p> <p>- Selection: Square wave (default setting)</p> <p>- Atrial Channel: Simulated P Wave</p> <p>- Width: 1 mS</p> <p>- Amplitude: 20 mVpeak</p> <p>- Active Outputs: Atrial channel (4mm banana jacks) only</p> <p>Square (SQU); triangle (TRI); haversine (HSN); sine square (SSQ); asymmetrical triangle (ISO) – fixed width: 2 mS rise time/13 mS fall time</p> <p>- Range: .05 mVpeak to 50.0 mVpeak</p> <p>- Accuracy: ± 5 % of setting</p> <p>- Resolution: 0.05 mVpeak (0.05 to 0.95 mVpeak); 0.50 mVpeak (1.0 to 49.5 mVpeak)</p> <p>- Range: 0.15 mS to 95.0 mS</p> <p>- Accuracy: ± 5 % of setting</p> <p>- Resolution: 0.05 mS (0.15 to 0.95 mS); 1 mS (1 to 19 mS); 5 mS (20 to 95 mS)</p> <p>Atrial channel (4mm banana jacks) only</p> <p>AAI (atrial pace and sense only)</p> <p>30 to 200 PPM</p> <p>500 Ω ± 1 %</p> <p>Paced refractory period; sensed refractory period</p> <p>20 to 500 mS</p> <p>± 5 % of reading (or ± 1 mS, whichever is greater)</p> <p>± 1 LSD</p> <p>3 digits</p> <p>- Selection: Square wave (default setting)</p> <p>- Ventricular Channel: Simulated R Wave</p> <p>- Width: 1 mS</p> <p>- Amplitude: 20 mVpeak</p>
--	---	---

	<p>- Active Outputs: Ventricular channel (4mm banana jacks) only</p> <p>Additional Waveform Selections: Square (SQU); triangle (TRI); haversine (HSN); sine square (SSQ); asymmetrical triangle (ISO) – fixed width: 2 mS rise time/13 mS fall time</p> <p>Amplitude:</p> <p>- Pacer Load Selection: 500 Ω</p> <p>-- Range: .05 mVpeak to 50.0 mVpeak</p> <p>-- Accuracy: $\pm 5\%$ of setting</p> <p>-- Resolution: 0.05 mVpeak (0.05 to 0.95 mVpeak); 0.50 mVpeak (1.0 to 49.5 mVpeak)</p> <p>- Default setting: 20 mVpeak</p> <p>Width:</p> <p>- Range: 0.15 mS to 300.0 mS</p> <p>- Accuracy: $\pm 5\%$ of setting</p> <p>- Resolution: 0.05 mS (0.15 to 0.95 mS); 1 mS (1 to 19 mS); 5 mS (20 to 95 mS); 25 mS (100 to 300 mS)</p> <p>- Default setting: 30 mS</p> <p>Intended Pacemaker Types: VVI</p> <p>Pacemaker Rates: 20 to 200 PPM</p> <p>DC Leakage Current</p> <p>Measurement Range: 00.1 μA to 99.9 μA</p> <p>Input Polarity: Positive and negative</p> <p>Resolution: 1 LSD (0.1 μA)</p> <p>Display Format: 3 digits</p> <p>Test Selections:</p> <p>- Static: Continuous (power OFF)</p> <p>- Dynamic: Gated (power ON)</p> <p>Test Load/Input Configurations:</p> <p>- Atrial+ and atrial-</p> <p>- Ventricular+ and ventricular-</p> <p>- Atrial+ and ventricular+</p> <p>Baseline/Test Selection: 500 Ω</p> <p>Dynamic Test Gating Algorithm: Measurement made 400 mS prior to the pacemaker pulse leading edge; 16 measurements averaged at a 4 mS rate for a total of 64 mS</p> <p>Specified Applied Pacemaker Rate: 80 PPM</p> <p>Current Drain Test</p> <p>DC Current Ranges: 0.100 mA to 0.999 mA; 1.00 mA to 9.99 mA; 10.0 mA to 99.9 mA</p> <p>Polarity: Positive or negative</p> <p>Indicator: + or – symbol</p> <p>Resolution: ± 1 LSD</p> <p>Display Format: 3 digits plus decimal point</p> <p>Accuracy: $\pm 5\%$ of reading $\pm 10 \mu$A</p> <p>Input DC Voltage:</p> <p>- Nominal: ± 9 V</p> <p>- Range: 5.0 V to 10.5 V</p> <p>- Input Protection: Short-circuit protection</p> <p>- Protection Type: Internal in-line fast-acting 1/2 A fuse</p> <p>Selectable Testloads: 200 Ω, 500 Ω, and 1000 Ω</p> <p>Battery Test Fixture: 9 V battery supply included, to facilitate connection of analyzer to recessed battery terminals within Medtronic 5388 and 5348 Temporary Pacemakers</p> <p>Test Loads</p> <p>Atrial Channel:</p> <p>- Selections: 200 Ω, 500 Ω, and 1000 Ω</p> <p>- Accuracy: $\pm 1\%$ of selection</p> <p>- Power Rating: 2 W</p> <p>Ventricular Channel:</p> <p>- Selections: 200 Ω, 500 Ω, and 1000 Ω</p> <p>- Accuracy: $\pm 1\%$ of selection</p> <p>- Power Rating: 2 W</p> <p>Tracking: Identical atrial and ventricular channel settings</p> <p>Input Defibrillation Protection:</p> <p>- Type: Internal spark gap</p> <p>- Episode Limit: 5 pulses @ 360 J (10 seconds min between discharges)</p> <p>- Life Limit: 250 pulses @ 360 J</p> <p>Long-Term Test</p>
--	--

	<p>Test Configuration:</p> <ul style="list-style-type: none"> - Transvenous Pacer: Atrial or ventricular channel only - Transcutaneous Pacer: Ventricular channel - Pulse Count Range: 999,999 (max) - Rate: 2 % to 20 % (default setting, 10 %) - Amplitude: 2 % to 20 % (default setting, 10 %) - Test Time (max): 999:59:59 (hhh:mm:ss) - Maximum Error Count: 200 - Test Termination: Manual; or upon max error count - Testloads: 200 Ω, 500 Ω, and 1000 Ω <p>Interactive Pacer ECG Simulation</p> <p>Simulation of demand, continuous, noncapture, and nonfunction patient</p> <p>Additional user</p> <p>-selectable parameters:</p> <ul style="list-style-type: none"> - NSR Heart Rate: Asystole and 20 to 250 BPM (1-BPM steps) - NSR PR Interval: 0.05 to 0.30 s (6 settings) <p>Pacemaker Capture/Threshold:</p> <ul style="list-style-type: none"> - Transcutaneous: 10 to 250 mA (10-mA steps) - Transvenous: 1 to 25 mA (1-mA steps) <p>General Specifications</p> <p>Serial Port:</p> <ul style="list-style-type: none"> - Type: RS232 - Connector Type: DB-9 (male) - Baud Rates: 2400, 9600, and 19200 - Data Control: Xon/Xoff <p>Power Requirements:</p> <ul style="list-style-type: none"> - External battery charger source/power supply - 100 to 240 Vac, 50/60 Hz operation - Auto power-off feature during battery operation <p>Battery Life:</p> <p>Environmental Specifications</p> <p>Temperature Range:</p> <ul style="list-style-type: none"> - Operating: 15 °C to 35 °C (59 °F to 95 °F) - Storage: 0 °C to 50 °C (32 °F to 122 °F) <p>Humidity Range:</p> <p>Dimensions</p> <p>Size:</p> <p>8" L X 4" W X 2" H (approx) (203 mm L X 101 mm W X 50 mm H)</p> <p>Weight:</p> <p>2 lb (approx) (0.90 kg)</p> <p>CE Mark</p> <p>User Safety:</p> <p>EMC:</p> <p>EN61326-1.1997</p> <p>Conforms to:</p> <p>UL STD 3101-1</p> <p>Certified to:</p> <p>CAN/USA STD C22.2 No.1010</p> <p>ETL Listed</p> <p>Device has received FDA 510(k) clearance (on file)</p>
--	--