

Maintenance Allocation Chart

4110-01-352-3653

Refrigerator, Mechanical, Blood Bank, Model FT2TRBLB

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

(1) GROUP NUMBER	(2) ASSEMBLY GROUP	(3) MAINTENANCE FUNCTION	(4) MAINTENANCE LEVEL	(5) TOOLS AND EQUIPMENT	(6) REMARKS
00	Refrigerator	Inspect	C		
		Test	C, O	19	
		Service	O	01, 02, 20, 29, 19	
		Repair	O	01, 02, 20, 29, 19	
		Replace	O, D	01, 02, 20, 29, 19	
		Overhaul	D	01, 02, 20, 29, 19	
		Rebuild	D	01, 02, 04, 20, 29, 19	
		01	Cabinet	Inspect	C
Test	C				
Service	C				
Repair	O			01, 02, 20, 29	
Replace	O			01, 02, 20, 29	
Overhaul	D			01, 02, 20, 29	
Rebuild	D			01, 02, 20, 29	
02	Light Switch			Inspect	C
		Test	C		
		Service	O	01, 02, 20, 29	
		Repair	O	01, 02, 20, 29	
		Replace	O	01, 02, 20, 29	
		Overhaul	O	01, 02, 20, 29	
		Rebuild	D	01, 02, 04, 20, 29	
		03	Light Bulb	Inspect	C
Test	C				
Service	C			01, 02, 20	
Replace	C			01, 02, 20	
04	Shelving	Inspect	C		
		Test	C		
		Service	C		
		Repair	O	01, 02	
		Replace	C	01, 02	
		Overhaul	O	01, 02	
		Rebuild	D	01, 02, 04	

Maintenance Allocation Chart

4110-01-352-3653

Refrigerator, Mechanical, Blood Bank, Model FT2TRBLB

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

(1) GROUP NUMBER	(2) ASSEMBLY GROUP	(3) MAINTENANCE FUNCTION	(4) MAINTENANCE LEVEL	(5) TOOLS AND EQUIPMENT	(6) REMARKS
05	Roller Drawers	Inspect	C		
		Test	C		
		Service	O	01, 02	
		Repair	O	01, 02	
		Replace	O	01, 02	
		Overhaul	O	01, 02	
		Rebuild	D	01, 02, 04	
06	Temperature Control	Inspect	C		
		Test	C		
		Service	O	01, 02, 20, 29, 19	
		Repair	O	01, 02, 20, 29, 19	
		Replace	O	01, 02, 20, 29, 19	
		Overhaul	O	01, 02, 20, 29, 19	
		Rebuild	D	01, 02, 04, 20, 29, 19	
07	Monitor Control Panel	Inspect	C		
		Test	O	01, 02	
		Service	O	01, 02, 20, 29, 19	
		Repair	O	01, 02, 20, 29, 19	
		Replace	O	01, 02, 20, 29, 19	
		Overhaul	O	01, 02, 20, 29, 19	
		Rebuild	D	01, 02, 04, 20, 29, 19	
08	Time Clock	Inspect	O		
		Test	O	01, 02	
		Service	O	01, 02, 29	
		Repair	O	01, 02, 29	
		Replace	O	01, 02, 29	
		Overhaul	O	01, 02, 29	
		Rebuild	D	01, 02, 04, 29	
09	Condensing Unit	Inspect	C	01, 02	
		Test	O	01, 02	
		Service	O	01, 02, 29	
		Repair	O	01, 02, 29	
		Replace	O	01, 02, 29	
		Overhaul	O	01, 02, 29	
		Rebuild	D	01, 02, 04, 29	

Maintenance Allocation Chart

4110-01-352-3653

Refrigerator, Mechanical, Blood Bank, Model FT2TRBLB

[C – Operator or Crew; O – Unit Maintenance; DS – Direct Support Maintenance; GS – General Support Maintenance; D – Depot Maintenance]

(1) GROUP NUMBER	(2) ASSEMBLY GROUP	(3) MAINTENANCE FUNCTION	(4) MAINTENANCE LEVEL	(5) TOOLS AND EQUIPMENT	(6) REMARKS
10	Blower Fan Assembly	Inspect	C		
		Test	O	01, 02,	
		Service	O	01, 02, 29	
		Repair	O	01, 02, 29	
		Replace	O	01, 02, 29	
		Overhaul	O	01, 02, 29	
		Rebuild	D	01, 02, 04, 29	
11	Defrost Heater	Inspect	O		
		Test	O	01, 02	
		Service	O	01, 02	
		Repair	O	01, 02, 29	
		Replace	O	01, 02, 29	
		Overhaul	O	01, 02, 29	
		Rebuild	D	01, 02, 04, 29	
12	Defrost Termination and Fan Switch	Inspect	O		
		Test	O	01, 02, 29, 19	
		Service	O	01, 02, 29, 19	
		Repair	O	01, 02, 29, 19	
		Replace	O	01, 02, 29, 19	
		Overhaul	O	01, 02, 29, 19	
		Rebuild	D	01, 02, 04, 29, 19	
13	Compressor	Inspect	O		
		Test	O	01, 02, 19	
		Service	D	01, 02, 04, 29, 19	
		Repair	D	01, 02, 04, 29, 19	
		Replace	D	01, 02, 04, 29, 19	
		Overhaul	D	01, 02, 04, 29, 19	
		Rebuild	D	01, 02, 04, 29, 19	