

DEPARTMENT OF THE ARMY SUPPLY BULLETIN

Army Medical Department Supply Information

Headquarters, Department of the Army, Washington, DC 20310-2300

20 April 2003

Effective until rescinded or superseded

DISTRIBUTION STATEMENT A: Approval for public release; distribution is unlimited

Printing of this information does not constitute endorsement of the product, manufacturer, or the use of any related item thereof by the U. S. Army.

REPORTING ERRORS AND RECOMMENDING IMPROVEMENTS

You can improve this publication. If you find any mistakes, or you know a way to improve it, please let us know. Mail your letter or DA Form 2028 (Recommended Changes to Publications and Blank Forms) directly to the U. S. Army Medical Materiel Agency, ATTN: MCMR-MMB-A, Fort Detrick, MD 21702-5001.

Table of Contents

Page

Overview of SB 8-75-S4	i
Chapter 1 - The AMEDD Logistics Assistance Program	1-1
Chapter 2 - Support Items of Equipment Common Deficiencies	2-1
Chapter 3 - Tactical Shelter Door Locking Device	3-1
Chapter 4 - M-400 Power Distribution Center PEU-155/E	4-1
Chapter 5 - Convention For Safe Containers (CSC)	5-1
Chapter 6 - Tactical Expandable Shelter Tiedown	6-1
Chapter 7 - Non-Medical ASIOE Repair Parts That Are Frequently Used For Repair	7-1
Chapter 8 - Assembly and Packing Instructions for Forward Surgical Team (FST) Operating Table	8-1
Chapter 9 - Tactical Shelter Scissors Jack Location for All Medical Materiel Sets (MMS)	9-1
Chapter 10 - The Acquisition Advice Codes (AAC) W and J	10-1
Chapter 11 - ISO and MILVAN Identification Markings	11-1
Chapter 12 - Operating Table Transport Restraint System, Model 2080M Surgical Tables	12-1
Chapter 13 - Complete Damage Inspection Required When Hospital Tactical Shelters and MILVANs are Moved	13-1
Chapter 14 - Air Conditioner Pulley/Bushing Problems	14-1
Chapter 15 - Unit Markings Placed on DEPMEDS Hospital Equipment	15-1
Chapter 16 - Inventory of DEPMEDS Equipment With TMs	16-1
Chapter 17 -	17-1
Chapter 18 -	18-1
Chapter 19 -	19-1
Glossary -	G-1
Index -	I-1

NOTICE

This Supply Bulletin is not superseded in its entirety.
Please refer to page i for supplemental information.

This Supply Bulletin is devoted entirely to
TOE Medical Units and Logistics Personnel.

SPECIAL NOTICE

Attention: TOE Medical Units and Logistics Personnel

Retain the 20 April 2001 and 20 April 2002 Editions of SB 8-75-S4.

The contents herein for 2003 are considered supplemental to the 2001 and 2002 editions. Or, consolidate all of 2002 and 2003 editions into the 2001 issue.

On a continuing bases, this SB will be updated either by including any revisions/changes in subsequent SBs, via our website <http://www.usamma.army.mil> or through the MMI message system. As events occur, we will keep all AIGs and SB addressees informed of the most current information available to us.

While not every page has been replaced or revised, there are changes that have occurred within the last year that are noteworthy to print. Please use this 2003 edition as a **supplemental resource - not a replacement or superseded edition**. The following list shows the status of the contents for 2003.

Chapter Number	CHAPTER TITLES	2003 STATUS: New (N) Replaced (R) Unchanged (U)
	Special Notice	R
	Overview	R
Chapter 1	The AMEDD Logistics Assistance Program	U
Chapter 2	Support Items of Equipment Common Deficiencies	U
Chapter 3	Tactical Shelter Door Locking Device	U
Chapter 4	M-400 Power Distribution Center PEU-155/E	U
Chapter 5	Convention For Safe Containers (CSC)	U
Chapter 6	Tactical Expandable Shelter Tiedown	U
Chapter 7	Non-Medical ASIOE Repair Parts That Are Frequently Used For Repair	U
Chapter 8	Assembly and Packing Instructions for Forward Surgical Team (FST) Operating Table	U
Chapter 9	Tactical Shelter Scissors Jack Location for All Medical Materiel Sets (MMS)	U
Chapter 10	The Acquisition Advice Codes (AAC) W and J	U
Chapter 11	ISO and MILVAN Identification Markings	U
Chapter 12	Operating Table Transport Restraint System, Model 2080M Surgical Tables	U
Chapter 13	Complete Damage Inspection Required When Hospital Tactical Shelters and MILVANs are Moved	U
Chapter 14	Air Conditioner Pulley/Bushing Problems	U
Chapter 15	Unit Markings Placed on DEPMEDS Hospital Equipment	U
Chapter 16	Inventory of DEPMEDS Equipment with TMs	U
Chapter 17	Maintenance and Repair of Tactical Shelter Lift and Support Jack Assemblies	U
Chapter 18	Problems With Temper Door Fasteners	U

(continued from page i)

Chapter Number	CHAPTER TITLES	2003 STATUS: New (N) Replaced (R) Unchanged (U)
Chapter 18	Problems With Temper Door Fasteners	U
Chapter 19	Hydrostatic Testing Of Compressed Gas Cylinders	U
Chapter 20	M10022A1 Dolly Set Status Update	U
Chapter 22	Management Of Non-Hospital Sets	U
Chapter 23	Augmentation Equipment And Its Associated Supplies	N
Glossary		R
Index		R

OVERVIEW

This Supply Bulletin is produced to assist the Logistics Personnel and the soldier in the field units in enhancing readiness with their equipment and within their Unit.

The feedback sheet below requests your proposals for improving the next edition of this Supply Bulletin. It also serves as a vehicle for submitting questions, problems and proposed solutions pertaining to non-medical equipment or ASIOE. The goal is to make future editions of this Supply Bulletin as informative and effective as possible.

CUSTOMER FEEDBACK For SB 8-75-S4 Dated 20 April 2003

Response From: _____

Telephone: _____ FAX: _____

E-mail: _____

To the USAMMA Website: <http://www.usamma.army.mil>

FEEDBACK (Please provide any constructive criticism about this edition):

COMMENTS (Please submit any questions, problems or solutions pertaining to non-medical equipment or ASIOE):

Send this sheet with comments to:

United States Army Medical Materiel Agency
 Attn: MCMR-MMR-M, Suite 100, Rm 104A
 1423 Sultan Dr.
 Fort Detrick MD 21702-5001
 DSN 343-4465 / Commercial 301-619-4465

OR

Contact our Customer Relationship Management Office (MCMR-MMB-R) from our website at <http://www.usamma.army.mil>, select Contact Us or also from the home page, under **Contact USAMMA**, click on the website and either send any additional feedback through the website facilities or call us at the phone number shown. We look forward to hearing from you.

CHAPTER 23. AUGMENTATION EQUIPMENT AND ITS ASSOCIATED SUPPLIES

a. In an attempt to provide the best medical care to the warfighter, the USAMMA has provided deploying units the most current technologies available on the commercial market. Some of these items are referred to as Augmentation items, defined as equipment that is not yet authorized, but has normally been documented by the Army Medical Department Center and School (USAMEDD C&S) as a requirement to increase the unit's ability to perform its medical mission. These items are currently being fielded to support current operations, and may not be fielded to all units as standard equipment.

b. These items will be issued to units via Letters of Authority (LOA) issued by the Major Command (MACOM). The LOA is the authorization document that allows the unit to retain the piece of equipment. Eventually the majority of these items will be documented to the units MTOE and the supporting supplies and accessories will be added to the unit assemblages.

c. These items are not documented in any unit assemblages, so it may be difficult to identify associated supplies and parts. The following information is provided to identify replacement items and to assist unit supply personnel in matching major end items to associated support equipment and parts.

1. OXYGEN CYLINDER REGULATOR "D" CYLINDER¹ – 6680-01-505-1303

¹ Allied Healthcare Products Model L270-220, all brass, oxygen regulator. Contains two DISS outlets at (40-60) psi with a constant flow of 100 lpm. Also includes a barb outlet for a cannula or mask with 12 flow settings ranging from 0 to 25 lpm. All three outlets can be used simultaneously; however, the maximum flowrate is 110 lpm divided among the three. Incorporates a CGA 870 connection for connection to a D-cylinder. Augments NSNs 6680-01-234-6789 and 6530-01-457-2239.

2. OXYGEN CYLINDER REGULATOR "H" CYLINDER² – 6680-01-505-1310

A review of the current hospital assemblages revealed a lack of 50-psi oxygen regulators required to interface oxygen H cylinders with the Impact 754M ventilator and the Narkomed M anesthesia machine. Additionally, it was determined that Level II and Level III medical treatment facilities did not possess enough of the D cylinder regulators. Based on this urgent need, the following regulators were identified for procurement.

² Allied Healthcare Products Model L280-220, all brass, oxygen regulator. Contains two DISS outlets at (40-60) psi with a constant flow of 100 lpm. Also includes a barb outlet for a cannula or mask with 12 flow settings ranging from 0 to 25 lpm. All three outlets can be used simultaneously; however, the maximum flowrate is 110 lpm divided among the three. Incorporates a CGA 540 connection for connection to an H-cylinder. Replaces NSN 6680-01-174-6276.

**3. PORTABLE OXYGEN GENERATION SYSTEM (POGS) –
6515-01-505-0203**

- Spare Parts Kit – System Level - 6515-01-504-1711
- Spare Parts Kit – MEDLOG BN/Depot Level – 6515-01-506-0009

**4. EXPEDITIONARY, DEPLOYABLE, OXYGEN CONCENTRATION
SYSTEM (EDOCS) – 6530-01-505-0526**

- Spare Parts Kit – System Level – Part No. 793018-001
- Spare Parts Kit – Organizational Level – Part No. 792923-001
- Spare Parts Kit – MEDLOG BN/Depot Level – Part No. 793017-00
- Patient Oxygen Distribution System (PODS) – 6530-01-447-4932
- Surgical Suite Oxygen Distribution System (SSODS) –
6530-01-447-4928
- Hospital Oxygen Backup System (HOBS) – 6530-01-447-8981

**5. COMPUTED RADIOGRAPHY READER (COMPANO) –
6525-01-505-7191**

- 8 X 10 Cassette – 6525-01-505-7232
- 10 X 12 Cassette – 6525-01-505-7230
- 14 X 17 Cassette – 6525-01-505-7228
- Un-interruptible Power Supply for Workstation – 6525-01-505-7221
- 15 inch LCD Color Monitor – 6525-01-505-7233
- 18 inch LCD Monitor – 6525-01-505-7234

**6. COMPUTED RADIOGRAPHY READER (OREX) -
6525-01-504-5002**

- CR 8 x 10 Plate and Cassette – Part No. 409081
- CR 10 x 12 Plate and Cassette Part No. 409080
- CR 14 x 17 Plate and Cassette Part No. 409079
- Toughbook 72 Military Laptop Part No. 409093
- Trip Lite UPS Model SU1000RT2U – Part No. PcUPS
- 650 MB Medical Grade CD-R in Jewel Box, 10/bx - Part No. 406179
- CD File, Write-On CD Bag with Hanger, 20per pack -
Part No. Boise S714403

**7. DENTAL DIGITAL IMAGER (DEXIS) –
6525-01-505-7780**

- Dexis Intra-oral sensor – Part No.DX-101
- Panasonic Toughbook 28 PIII 800 Laptop – Part No. CF-28PBJGZPM

**8. ULTRASOUND, HAND HELD (SONOSITE 180 Plus) –
6525-01-503-7182/FST and 6525-01-503-7170/HOSP**

- Microcurved Array C15 4/2 - 6515-01-502-3728
- Intrecavatory Array ICT 7-4 - 6515-01-502-3734
- Curved Array C60 5-2 - 6515-01-502-3739
- Site Charge Daul Battery Charger - 6130-01-490-2678
- Latex Probe Cover - 6515-01-504-6163

GLOSSARY FOR SB 8-75-S4, 2003

Acronym	Explanation
AAC	Acquisition Advice Code
AAC-L	Acquisition Advice Code Local Purchase non-Stocked Items
ACSLOG	Assistant Chief of Staff for Logistics
AMEDD	Army Medical Department
AMEDDC&S	Army Medical Department Center & School
AR	Army Regulation
ARNG	Army National Guard
ASIOE	Associated Support Items of Equipment
ATRRS	Army Training Requirements and Resources System
BII	Basic Issue Items
CAGE	Contractor/Commercial and Government Entity
CARC	Chemical Agent Resistance Compound
CATC	Combined Arms Training Center
CONUS	Continental United States
CSC	Container Safety Certifications, Convention for Safe Containers
CSDP	Command Supply Discipline Program
CSH	Combat Support Hospital
DA	Department of the Army
DCMC	Defense Contract Management Command
DD	Defense Department
DEPMEDS	Deployable Medical Systems
DES	Dental Equipment Sets
DISE	Distribution Illumination Systems Electrical
DOD	Department of Defense
DOL	Directorate of Logistics
DOT	United States Department of Transportation
DPSC	Defense Personnel Support Center
DRMO	Defense Reutilization & Marking Office
DRMS	Defense Reutilization & Marking Service
DCSOPS	Deputy Chief of Staff for Operations and Plans
DTS	Defense Transportation System
DX	Direct Exchange
EAC	Echelon Above Corps
EC	Essential Characteristics
ECU	Environmental Control Unit
ERC	Equipment Readiness Code
ESCO	Engineered Systems Company
FDA	Food and Drug Administration
FEA	Finite Element Analysis
FEDLOG	Federal Logistics Record
FMC	Fully Mission Capable
FST	Forward Surgical Team

(continued) GLOSSARY FOR SB 8-75-S4, 2003

Acronym	Explanation
GPM	Ground Precautionary Message
GS	General Support
HAZ-MAT	Hazardous Material
IAW	In Accordance With
IMM	Integrated Materiel Manager
ISO	International Standardization Organization
IV	Intravenous
JRTC	Joint Readiness Training Center
KRC	Keweenaw Research Center
kw	Kilowatts
LAP	Logistics Assistance Program
LAR	Logistical Assistance Representative
LAV	Logistics Assistance Visit
LIN	Line Item Number
LO	Lubrication Order
LOA	Letter of Authority
MACOM	Major Command
MBU	Modern Burner Unit
MEDLOG	Medical Logistics System
MESs	Medical Equipment Sets
MF2K	Medical Force 2000
MILVAN	Military-Owned De-mountable Container
MMS	Medical Materiel Sets
MOTSU	Military Ocean Terminal Sunny Point
MRI	Medical Reengineering Initiative
MTT	Mobile Training Team
MWO	Modification Work Order
NRC	Nonreusable, Nonrefillable Containers
NSN	National Stock Number
NTC	National Training Center
OCONUS	Outside Continental United States
OMA	Operation and Maintenance, Army
OPA	Other Procurement, Army
OR	Operating Room
PMCS	Preventive Maintenance Checks and Services
POC	Point of Contact
PUP	Polyurethane Paint
QDR	Quality Deficiency Report

(continued) GLOSSARY FOR SB 8-75-S4, 2003

Acronym	Explanation
RPSTL	Repair Parts and Special Tools List
RTCH	Rough Terrain Cargo Handler
RTS-MED	Reserve Training Site Medical
SMART	Special Management Augmentation Reaction Teams
SOUN	Safety of Use Message
TACOM	Tank Automotive Command
TARDEC	TACOM Research and Development Center
TDP	TACOM Drawing Package
TDY	Temporary Duty
TEMPER	Tent, Extendable, Modular, Personnel
TM	Technical Manual
TMDE	Test, Measurement and Diagnostic Equipment
TT	Total Time
UA	Unit Assemblage
UAL	Unit Assemblage Listing
USACHPPM	U.S. Army Center for Health Promotion and Preventive Medicine
USAH	United States Army Hospital
USAH/U	United States Army Hospital / Unit
USAMMA	United States Army Medical Materiel Agency
USAREUR	United States Army, Europe
USAMEDD C&S	United States Army Medical Department Center & School
USP	United States Pharmacopoeia
VES	Veterinary Equipment Sets
VHS	Vertical Helix Scan
VOC	Volatile Organic Concentrator

INDEX FOR SB 8-75-S4, 2003


<u>SUBJECT</u>	<u>SB 8-75-</u>	<u>PAGE</u>
Acquisition Advice Code (AAC) Of W & J and How They are Used	S4	10-1
Air Conditioner Pulley/Bushing Problems	S4	14-1
Air Conditioner, AE32C39	S4	2-18
Assembly Instructions for OR Table (NSN 6530-01-321-5592)	S4	8-2
Augmentation Equipment And Its Associated Supplies	S4	23-1
CARC, Lessons Learned	S4	15-1
Compliance Requirements for the International Safety Standards to Certify Containers	S4	5-1
Component Descriptions for OR Table	S4	8-10
Components of OR Table (NSN 6530-01-321-5592)	S4	8-1
Container, General Cargo, Military-Owned Demountable Container (MILVAN)	S4	2-6
Containers Require the Proper Markings Prior to Recertification	S4	11-1
Course Schedule for FY 2002 Intermodal Dry Cargo Container Reinspection	S4	5-2
DA Form 2404 (Sample Form)	S4	Appendix B
DD Form 2282 (Sample Form)	S4	Appendix C
DEPMEDS Equipment is Being Painted with Other than CARC Paint	S4	15-1
Description of Table, Operating, Field (NSN 6530-01-321-5592)	S4	8-1
Disassembly and Repacking Instructions for Field Operating Table (NSN 6530-01-321-5592)	S4	8-18
DOD Form DD 2282, Transmittal Request for New	S4	5-4
Dolly Set, M-1022 (NSN 2330-01-167-7262)	S4	2-7
Dolly Set, M-1022A1 (NSN 2330-01-378-9997)	S4	2-7
Fan Pulleys on Air Conditioner Units Present Problems	S4	14-1
Fielding Status of Modern Burner Unit to the TOE Hospital Units	S4	21-2
Forward Surgical Team Operating Table Assembly Instructions	S4	8-1
Frequently Used Repair Parts for the Non-Medical Equipment	S4	7-1
Generator, 100 Kilowatts (kw)	S4	2-17
Health/Safety/Environmental Hazards of CARC	S4	15-2
Improvement of 10/20 Maintenance Standards at Unit Level	S4	2-1
Information on the History and Current Issues with the M1022A1 Dolly Set	S4	20-1
Installation Maintenance, Lubrication and Repair of Tactical ISO Shelter Lift and Support Jack Assemblies	S4	17-1
Inventory of DEPMEDS Equipment is Difficult when Using TMs	S4	16-1
ISO Jack(s) Requirement(s)	S4	9-1
Jack Assembly Hinged (Figs. 17-1 through 17-3)	S4	17-2
Lessons Learned about CARC	S4	15-1
Logistics Assistance Program (LAP), The	S4	1-1
M1022A1 Dolly Set, Cracks and Location of Cracks (Fig. 20-1)	S4	20-2
M1022A1 Dolly Set, Information on the History and Current Issues of	S4	20-1
Modern Burner Unit (Fig. 21-1)	S4	21-1
Modern Burner Unit, Fielding Status to the TOE Hospital Units (Table 21-1)	S4	21-2

<u>SUBJECT</u>	<u>SB 8-75-</u>	<u>PAGE</u>
Non-Hospital Set Management, Update on Current Issues	S4	22-1
Non-Hospital Sets, Potential Change in USAMRMC Fielding Policy	S4	22-1
Nonmedical Associated Support Items of Equipment (ASIOE) Repair Parts that are Frequently Used to Repair the ASIOE	S4	7-1
Non-Medical Equipment Inspected During Field Visits	S4	2-1
Operating Table Transport System Model 2080M Surgical Tables	S4	12-1
OR Table (NSN 6530-01-321-5592)	S4	8-15
Panel, Power Distribution, M100	S4	2-20
Panel, Power Distribution, M40	S4	2-19
Panel, Power Distribution, M46	S4	2-20
Panel, Power Distribution, M4-600, PEU-155/E	S4	2-19
Parts Listing for Non-Medical Equipment (ASIOE)	S4	7-1
Potential Change in U.S. Army Medical Research and Materiel Command (USAMRMC) Fielding Policy for Non-Hospital Sets	S4	22-1
Proper Care and Maintenance of the Slide Fasteners on the Medical TEMPER	S4	18-1
Proper Placement of Scissors Jack(s) to Support the ISO Shelter	S4	9-1
Proper Tiedown Procedures	S4	6-2
Refrigerated MILVAN	S4	2-22
Safety Inspection Requirements for CSC	S4	5-1
Security and Safety of the Tactical Shelter	S4	3-1
Shelters, Tactical Expandable	S4	2-2
Storage Container (NSN 6530-01-432-3173) for OR Table (NSN 6530-01-432-3173)	S4	8-18
Tactical Shelter Door Locking Device - Introduction	S4	3-1
Tactical Shelter Scissors Jack Location for Added Support	S4	9-1
Tactical Shelter Security and Safety	S4	3-1
Tactical Shelters & MILVANs Require Inspection when Moved	S4	13-1
Tactical Shelters and MILVANs Required Safety Inspection	S4	5-1
Tactical Shelters Damaged from Improper Tie-down	S4	6-1
TEMPER 6x20 Surgical Work Area (Fig. 18-1)	S4	18-1
TEMPER Slide Fasteners, Proper Care of	S4	18-1
Tent, Expandable, Modular, Personnel (TEMPER)	S4	2-21
Testing and Proper Storage of Medical Oxygen Cylinder	S4	19-1
The Modern Burner Unit is the Replacement for the M2 Gasoline Burner Currently Used in All Field Feeding Systems to Include DEPMEDS	S4	21-1
TM43-0139, Excerpt About CARC Hazards (Fig. 15-1)	S4	15-2
Transmittal Request for New DOD Form DD 2282	S4	5-4
Update on Current Issues with Non-Hospital Set Management	S4	22-1
Updated Status on the MWO and Retrofit Kits	S4	20-2
Use of Additional Support to the Tactical ISO Shelter	S4	9-1
Wiring Chaffing Stripe Missing Inside The M-400 Case	S4	4-1

By Order of the Secretary of the Army:

ERIC K. SHINSEKI
General, United States Army
Chief of Staff

Official:


JOEL B. HUDSON
Administrative Assistant to the
Secretary of the Army

Distribution:

To be distributed in accordance with initial distribution number (IDN) 340016, requirements for SB 8-75 Series.